

Neo-Flex^{ACTIVE}TM

Flexed Suspension Sleeve




You Customize the Length!

The Neo-Flex Active is a durable, cost-effective prosthetic suspension sleeve. It is constructed of a strong, latex-free neoprene with built-in knee pad for added cushioning and wear. Its flexed shape and Nylon knee center add softness and flexibility to enhance comfort and mobility. The Neo-Flex Active may be trimmed at proximal and distal ends to customize length if needed.

Trimming Instructions:

If length adjustment of proximal and/or distal end of sleeve is desired, choose any spot BETWEEN the reinforcement tacks located along the posterior seam of the sleeve and mark cutting line with chalk or marker. Place sleeve on its side and cut from posterior to anterior along your cutting line using sharp scissors.



Indications:

For use as a prosthetic suspension aid.

Warning:

Do Not use over broken or irritated skin. Discontinue use of sleeve if any skin irritation, blisters or rash develop with use. Do not wear while sleeping or use if you have circulatory problems.

Care Instructions:

Wash sleeve with mild detergent and cool water and lay flat to dry. Do not use excessive heat, chlorine bleach, or solvents.

Warranty:

Defects in manufacturing or materials used will be covered for 90 days from the date of purchase. Defective item must be returned clean and dry and will either be replaced or repaired.

Neo-Flex^{ACTIVE}TM Suspension Sleeve Sizing

For proper sizing, take thigh circumference measurement 5" proximal to top edge of the patella and refer to sizing chart.

Neo-Flex Active is available in a two-tone Gray/Black.

Limb Measurement	Sleeve Size	Ref no.
10.5" -12.5"	0	CA700-FPA-0
11.5" -13"	1	CA700-FPA-1
12.5" -14"	1.5	CA700-FPA-15
14" -16"	2	CA700-FPA-2
15" -17"	2.5	CA700-FPA-25
16" -18"	3	CA700-FPA-3
18.5" -20	4	CA700-FPA-4
21" -23"	5	CA700-FPA-5